

## Assembly Wrench for UW 50 and LS 20 - Bandelin Sonopuls AH 50

Item no. BND-3820

Manufacturer Bandelin

AH 50 assembly tool for mounting and removing the UW 50 ultrasonic transducer and the LS 20 lab rod. It allows the components to be screwed together safely and without damage on the SONOPULS homogenizer.

### TECHNICAL DATA

|                        |          |
|------------------------|----------|
| Ausführungsklasse      | TEC      |
| Country of Manufacture | Germany  |
| HS-Code                | 84799080 |
| Weight                 | 0.184 kg |



### DESCRIPTION

The AH 50 is the assembly tool for fastening and releasing the UW 50 ultrasonic transducer and the LS 20 lab rod. It allows the components to be screwed together in a defined way without damaging the sensitive threads and sound-transmitting faces.

### BENEFITS

- Fits the UW 50 ultrasonic transducer and the LS 20 lab rod
- Safe assembly without damage to threads or sound faces
- Eases the changing of sonotrodes and horns

### TECHNICAL DATA

|                 |              |
|-----------------|--------------|
| Order no.       | 3820         |
| Compatible with | UW 50, LS 20 |
| Weight          | 0.184 kg     |

### APPLICATIONS

The AH 50 is used in the lab when setting up and reconfiguring the SONOPULS homogenizer with the UW 50 and LS 20, for example when changing sonotrodes or stepped horns.