

Glass Beaker 600 ml for Indirect Sonication - Bandelin Sonorex SD 05

Item no. BND-575

Manufacturer Bandelin

Glass insert beaker of 600 ml for hanging into a perforated lid or insert strip. It enables indirect sonication of aggressive or sensitive media without loading the oscillating tank.

TECHNICAL DATA

Ausführungsklasse	TEC
Country of Manufacture	Germany
HS-Code	84799080
Weight	0.757 kg



DESCRIPTION

The SD 05 is a glass insert beaker with 600 ml capacity for indirect sonication. Using the GR 04 ring, the beaker is hung into a perforated lid or insert strip so that aggressive cleaners, abrasive soils or sample liquids can be sonicated without anything reaching the tank floor.

BENEFITS

- 600 ml glass beaker for indirect sonication
- Hangs into a perforated lid or insert strip with the GR 04 ring
- Several beakers with different liquids can be used at once
- Protects the oscillating tank from aggressive media

TECHNICAL DATA

Order no.	575
Working volume	0.6 L
Inner dimensions (L×W×D)	76 × 76 × 150 mm
Outer dimensions (L×W×H)	80 × 80 × 150 mm
Weight	0.757 kg

APPLICATIONS

For example, dental cements can be treated with STAMMOPUR Z in one beaker and dental instruments with STAMMOPUR DR 8 in a second beaker at the same time. To transmit the ultrasound the surface tension of the bath water must be reduced, so a cleaning solution such as TICKOPUR R 33 at 3 percent is used as coupling liquid.

Compatible with the DE 08 perforated lid.