

Sieve Basket with Fine Mesh for Insert Vessels - Bandelin Sonorex PD 04

Item no. BND-126

Manufacturer Bandelin

Fine-mesh sieve basket for use in insert vessels or for cleaning very fine parts directly in the ultrasonic bath. The close mesh holds the smallest parts securely while letting the cleaning fluid pass freely.

TECHNICAL DATA

Ausführungsklasse	TEC
Country of Manufacture	Germany
HS-Code	84799080
Weight	0.126 kg



DESCRIPTION

The PD 04 sieve basket securely holds very fine and small parts during ultrasonic cleaning. The close mesh lets the cleaning fluid and cavitation pass through unhindered, so even the smallest components are cleaned thoroughly and are not lost through the basket.

BENEFITS

- Fine mesh holds the smallest parts securely
- Cleaning fluid circulates freely through the sieve
- For use in insert vessels or directly in the bath
- Compact form for targeted small-part cleaning

TECHNICAL DATA

Order no.	126
Inner dimensions (L×W×D)	60 × 60 × 45 mm
Outer dimensions (L×W×H)	75 × 75 × 50 mm
Weight	0.126 kg

APPLICATIONS

The sieve basket is suited to cleaning very fine parts such as screws, seals, small components and

jewellery pieces, either directly in the ultrasonic bath or placed inside an insert vessel for indirect sonication.

The PD 04 fits the insert vessels KB 04, SD 04, SD 05, SD 06, SD 09, EB 05 and PD 06.