

Perforated Lid for Insert Beakers, RK/DT 100 to 103 - Bandelin Sonorex DE 100

Item no. BND-3017

Manufacturer Bandelin

Stainless steel perforated lid that holds insert beakers for the 100 to 103 series ultrasonic baths. It enables indirect sonication of chemically aggressive or abrasive media in the beaker.

TECHNICAL DATA

Ausführungsklasse	TEC
Country of Manufacture	Germany
HS-Code	84799080
Weight	0.516 kg



DESCRIPTION

The DE 100 perforated lid holds insert beakers above the oscillating tank and so enables indirect sonication. Chemically aggressive liquids, abrasive contamination or sample fluids are treated inside the beaker without loading the tank floor, since nothing may be placed directly on the bottom.

BENEFITS

- Indirect sonication of delicate or aggressive media in the insert beaker
- Protects the tank floor, as there is no direct floor contact
- Several beakers with different liquids can be sonicated at once
- Corrosion-resistant stainless steel, easy to clean

TECHNICAL DATA

Order no.	3017
Outer dimensions (L×W×H)	264 × 132 × 11 mm
Tank material	Stainless steel
Weight	0.516 kg

APPLICATIONS

Typical uses are removing dental cements with STAMMOPUR Z, chemically disinfecting dental

instruments with STAMMOPUR DR 8, or treating samples with TICKOPUR J 80 U, each inside the insert beaker and without contaminating the bath fluid.

Fits the ultrasonic baths RK 100/H, RK 102 H, RK 103 H, DT 100/H, DT 102 H/H-RC, DT 103 H and DL 102 H.